

## Accepted Manuscript

Title: Capillary thermal desorption unit for near real-time analysis of VOCs at sub-trace levels. Application to the analysis of environmental air contamination and breath samples

Authors: Mónica Alonso, Mar Castellanos, José Martín, Juan M. Sanchez

PII: S1570-0232(09)00202-5  
DOI: doi:10.1016/j.jchromb.2009.03.022  
Reference: CHROMB 16271

To appear in: *Journal of Chromatography B*

Received date: 15-12-2008  
Revised date: 17-3-2009  
Accepted date: 17-3-2009

Please cite this article as: M. Alonso, M. Castellanos, J. Martín, J.M. Sanchez, Capillary thermal desorption unit for near real-time analysis of VOCs at sub-trace levels. Application to the analysis of environmental air contamination and breath samples, *Journal of Chromatography B* (2008), doi:10.1016/j.jchromb.2009.03.022

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



1 **Capillary thermal desorption unit for near real-time analysis**  
2  
3  
4 **of VOCs at sub-trace levels. Application to the analysis of**  
5  
6  
7 **environmental air contamination and breath samples**  
8  
9

10  
11  
12  
13 **Mónica Alonso<sup>1</sup>, Mar Castellanos<sup>2,3</sup>, José Martín<sup>1</sup> and Juan M. Sanchez<sup>1,3,\*</sup>**  
14  
15

16  
17  
18  
19 <sup>1</sup> Department of Chemistry, University of Girona, Girona (Spain)  
20

21  
22 <sup>2</sup> Department of Neurology, Dr Josep Trueta University Hospital, Girona (Spain)  
23

24  
25 <sup>3</sup> Girona Biomedical Research Institute (IdiBGi), Girona (Spain)  
26  
27  
28  
29  
30  
31  
32

33 \* Corresponding author: [juanma.sanchez@udg.edu](mailto:juanma.sanchez@udg.edu), Phone (+34) 972418276, Fax:  
34  
35 (+34) 972418150  
36  
37  
38  
39  
40

41 **Keywords: thermal desorption, breath, smoking, biomarker**  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

## ***Abstract***

1  
2  
3 A capillary microtrap thermal desorption module is developed for near real-time  
4 analysis of volatile organic compounds (VOCs) at sub-ppbv levels in air samples. The  
5 device allows the direct injection of the thermally desorbed VOCs into a  
6 chromatographic column. It does not use a second cryotrap to focalize the adsorbed  
7 compounds before entering the separation column so reducing the formation of artifacts.  
8 The connection of the microtrap to a GC-MS allows the quantitative determination of  
9 VOCs in less than 40 minutes with detection limits of between 5-10 pptv (25°C and 760  
10 mm Hg), which correspond to 19-43 ng·m<sup>-3</sup>, using sampling volumes of 775 cm<sup>3</sup>. The  
11 microtrap is applied to the analysis of environmental air contamination in different  
12 laboratories of our faculty. The results obtained indicate that most volatile compounds  
13 are easily diffused through the air and that they also may contaminate the surrounding  
14 areas when the habitual safety precautions (e.g., working under fume hoods) are used  
15 during the manipulation of solvents. The application of the microtrap to the analysis of  
16 VOCs in breath samples suggest that 2,5-dimethylfuran may be a strong indicator of a  
17 person's smoking status.  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42

## ***Introduction***

43  
44 The concept that blood, urine, and other body fluids and tissues can be collected and  
45 analyzed to yield information for the diagnosis of disease states or to monitor disease  
46 progression and/or therapy is an essential concept in modern medicine. Many important  
47 developments in medical monitoring technologies and diagnostic methods have focused  
48 on blood and urine analysis for clinical diagnosis. Diagnoses based on breath analysis  
49 are much less developed and not yet widely used in clinical practice [1], in spite of the  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 fact that breath measurement has enormous potential, especially given its inherent  
2 safety. The only requirement for the collection of a breath sample is that the subject  
3 must be breathing. However, despite its great promise and the fact that there have been  
4 hundreds of publications on breath analysis there are only a handful of tests that are  
5 used clinically and a few others that are used in research [2].  
6  
7  
8  
9  
10

11  
12  
13  
14  
15 Breath testing challenges scientists in two ways. First, it is technically very difficult to  
16 analyze breath volatile organic compounds (VOCs) that are present at picomolar  
17 concentrations. Second, the source and biochemical significance of most breath VOCs  
18 are still unknown [3,4]. It is important to develop methods capable of trace  
19 determination and so a preconcentration technique is required before the introduction of  
20 the collected sample into a gas chromatographic column. To meet the need for clinical  
21 applications, a relatively inexpensive, portable instrument capable of providing  
22 nonintrusive, real-time, sensitive, and accurate analysis of breath gases for medical  
23 diagnosis is highly desirable.  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39

40 In current public health, smoking is one of the most important single risk factors and is  
41 related to several pathological conditions, mainly pulmonary diseases. However, so far  
42 we have not identified any valid biomarker that serve as proxies for tobacco exposure  
43 and their possible association with and predictive capacity of tobacco-related diseases.  
44  
45  
46  
47  
48  
49 This highlights the importance and urgency of conducting research in this area [5].  
50  
51  
52  
53  
54

55 A variety of biomarkers to evaluate environmental tobacco smoke (ETS) exposure have  
56 been proposed. Different chemicals present in tobacco smoke, such as carbon monoxide  
57 and cyanide (the latter metabolized in the body to thiocyanate), can be analyzed in  
58  
59  
60  
61  
62  
63  
64  
65

1 biological fluids. However, these chemicals are not specific markers of ETS exposure as  
2 there are significant sources of carbon monoxide and cyanide, including the body's own  
3 metabolism, other than ETS. The measurement of nicotine metabolites in biological  
4 fluids (mainly blood and urine) has attracted more attention from the scientific  
5 community. The determination of cotinine, the major proximate metabolite of nicotine,  
6 has been widely used by scientist to evaluate ETS exposure given that this substance  
7 reflects exposure to nicotine, which is almost specific to tobacco [6]. Nicotine-derived  
8 nitrosoamines such as 4-(methylnitrosoamino)-1-(3-pyridyl)-1-butanone (NNK) have  
9 also been proposed [7]. This nitrosoamine is metabolized to a butanol metabolite (4-  
10 (methylnitrosoamino)-1-(3-pyridyl)-1-butanol, NNAL), and its glucuronide (NNAL-  
11 GLUC). From the different metabolites evaluated, cotinine appears to be the most  
12 specific and most sensitive biomarker for exposure to nicotine from ETS. A limitation  
13 of using cotinine, however, is that it indicates ongoing exposure but not long-term  
14 exposure to ETS. A further limitation to the use of nicotine metabolites in general is that  
15 their determination is complex, requiring the taking of urine or blood samples. The  
16 possibility of using breath tests to evaluate exposure to ETS is very promising given  
17 that it is a procedure that is safe and easy to perform. Preliminary studies have  
18 suggested that some VOCs detected in breath samples may be biomarkers (e.g., 1,3-  
19 butadiene, benzene, 2,5-dimethylfuran) [8-16].

20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49 The aim of this study is to develop a capillary microtrap thermal desorption module  
50 without the use of a secondary cryogenic trap for the focusing of the retained  
51 compounds before entering the column in the GC system. This in-house device has been  
52 specifically designed for analysis of VOCs at ppbv-pptv levels in breath samples. The  
53 system allows fast analysis (less than 40 minutes) of air sample volumes below 1 L.  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 This will help to monitor near real-time environmental contamination and it will be used  
2 for the study and evaluation of smoking biomarkers in breath. The concentration device  
3  
4 used here is a variation of one previously described [14,17].  
5  
6  
7  
8  
9  
10

## 11 ***Experimental***

### 12 13 14 15 16 17 18 **Apparatus**

19  
20  
21  
22  
23 The experimental system used for these studies is shown in Figure 1. There is currently  
24 no commercial instrument able to analyze sub-ppmv VOCs in breath samples on-line. It  
25 is therefore interesting to develop an in-house concentrator/injector device specifically  
26 designed for this type of analysis. Such a device should be able to perform two main  
27 functions: firstly, to retain target VOCs quantitatively from samples and, secondly, to  
28 achieve a fast and quantitative desorption of the adsorbed analytes to be sent to a GC  
29 column as a small injection band, preventing the formation of broad peaks with tails in  
30 the chromatograms.  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44

45 A previous study [17] describes in detail the design of a similar microtrap configuration  
46 as the one used in the present study, and the heating process and control circuitry have  
47 also been described. As there is no single sorbent material that meets the full ideal  
48 criteria for concentrating all VOCs from an air matrix, a multi-bed trap must be used for  
49 this purpose [17,18]. In this study, we used a three-bed trap for general breath  
50 monitoring as it has been found to be effective in analyzing common VOCs present in  
51 breath samples [14]. Quantities of  $2.5 \pm 0.2$  mg of Carboxen 1000 and Carboxen X  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 (Supelco, Bellefonte, PA, USA) and  $5.5 \pm 0.2$  mg of Carbopack B (Supelco) were  
2 sequentially introduced in an 80 mm long, 1.35 mm ID Ni/Co alloy tube (Accu-Tube  
3 Corp., Englewood, CO, USA). The packed tube was preconditioned by heating at  $300^\circ\text{C}$   
4 for 4 h with a continuous nitrogen flow ( $200 \text{ mL} \cdot \text{min}^{-1}$ ). Sampling and desorption flow  
5 were in opposite directions to avoid the least volatile compounds entering into contact  
6 with the strongest adsorbent (Carboxen 1000). Under the conditions employed,  
7 desorption was complete and no memory effects were observed in subsequent heating  
8 cycles.  
9

10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21 A vacuum pump (Air Cadet Vacuum Station, Barnant Co., Barrington, IL, USA) was  
22 used to pull sample gas from left to right through the trap tube shown in Figure 1 at a  
23 sample flow rate of  $31 \text{ cm}^3 \cdot \text{min}^{-1}$ . A three-way valve V1 (01380-05, Cole Parmer,  
24 Vernon Hills, IL, USA) and two-way valves V2 and V3 (LFVA1230113H, Lee Co.,  
25 Westbrook, CT, USA) were used for flow control of the sample and the carrier gas. For  
26 sample collection, V2 and V3 were opened and V1 was set to the vacuum pump  
27 direction. After sample collection was complete, valve V2 was closed for 10 s to let all  
28 the sample gas remaining in the lines go through the trap. Valve V3 was then closed and  
29 valve V1 was switched to allow the carrier gas to move through the trap line. The use of  
30 a bypass carrier gas line (valve V3 line) was necessary during sample collection as  
31 otherwise part of the sample moves directly to the separation columns due to the strong  
32 vacuum pump of the MS detector. In the case of breath samples, this problem is  
33 aggravated due to the presence of substantial amounts of water vapor. A previous study  
34 [14] shows that this configuration helps to maintain the pneumatic conditions in the  
35 separation column during the sampling process and prevent thermal decomposition on  
36 the surface of the carbon-based sorbents of the most reactive compounds in the samples.  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1  
2 Monitoring of the desorption temperature was performed with a J-type thermocouple  
3  
4 using 0.127 mm ID wire (36 AWG, Omega Engineering, Stamford, CT, USA) attached  
5  
6 to the outer wall of the tube. In order to obtain very rapid heating with a minimal  
7  
8 overshoot, a current pulse was applied to the trap for 0.9 s increasing the temperature  
9  
10 from 25°C to the desired desorption temperature (~300°C). The use of a second low  
11  
12 current pulse to maintain the trap temperature, as described in previous studies [14,17],  
13  
14 was not necessary as all the adsorbed compounds are quantitatively desorbed with the  
15  
16 first pulse.  
17  
18  
19  
20  
21  
22  
23  
24  
25

## 26 **Materials and methods**

27  
28  
29  
30  
31 All analytes were reagent grade or better. Standards were prepared in 10 L Tedlar gas  
32  
33 sampling bags (SKC Inc., Eighty Four, PA, USA) by injecting 1 to 2  $\mu\text{L}$  of individual  
34  
35 components and diluting them with dry air. Vapor concentrations were in the range of  
36  
37 8–46 ppmv.  
38  
39  
40  
41  
42

43  
44 Component separation was achieved by the use of a 30 m length of nonpolar 5% phenyl  
45  
46 95% dimethylpolysiloxane column with 0.25 mm ID and 0.25  $\mu\text{m}$  film thickness (DB-5,  
47  
48 J&W Scientific, Folsom, CA, USA for FID applications, and ZB-5ms, Zebron,  
49  
50 Phenomenex, Torrance, CA, USA for MS applications). For the evaluation of the  
51  
52 developed on-line microtrap, a GC 8000 series (Fisons Instruments, Milano, Italy) with  
53  
54 an FID at 250°C was used. For the analysis of the breath and air samples, a Focus GC  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65



1 (Thermo Scientific, Waltham, MA, USA) with a mass spectrometer detector (DSQ II,  
2 Thermo Scientific) was used.  
3

4  
5  
6  
7 Helium carrier gas was used after purification for water vapor, hydrocarbons, and  
8  
9 oxygen. Instrument and data control was provided by a 12-bit A/D board (PCI-1710HG,  
10  
11 Advantech, Taiwan). A sampling rate of 10 Hz was selected.  
12  
13

### 14 15 16 17 18 19 **Sample Collection** 20

21  
22  
23  
24 Breath samples from students and staff of the University of Girona were collected in 1 L  
25  
26 Tedlar gas-sampling bags and analyzed immediately after collection. For each sample,  
27  
28 775 cm<sup>3</sup> of breath was pulled through the multibed sorption microtrap at a rate of 31  
29  
30 cm<sup>3</sup>·min<sup>-1</sup>. It is known that significant changes in breath molecule concentrations take  
31  
32 place as a result of variations in ventilation parameters during sampling, and  
33  
34 hyperventilation may occur when a breath sample is collected from spontaneously  
35  
36 breathing subjects [19]. However, this study was not aimed at determining the amounts  
37  
38 of compounds from metabolic sources but rather at developing a reliable and  
39  
40 quantitative method with low detection limits that allows the determination of adequate  
41  
42 breath biomarkers to establish smoking status. For this reason, spontaneously breathing  
43  
44 subjects were used. During exhalation, the first 5 s was not collected in the sampling  
45  
46 bag in order to reduce the contamination by dead-space air.  
47  
48  
49  
50  
51  
52

53  
54  
55 For the analysis of environmental air, 1 L of sample was taken with a 1 L gas syringe  
56  
57 (SGE JUMBO syringe, SGE Europe, UK), the sample was introduced in a 1 L Tedlar  
58  
59  
60  
61  
62  
63

1 bag and analyzed with the system developed. The same sampling volume as in the  
2 breath analysis was used unless otherwise specified. Smaller sample volumes were used  
3  
4 in situations of large contamination, where sampling time was decreased to avoid the  
5  
6 breakthrough of sorbents.  
7  
8  
9

## 10 11 12 13 14 ***Results and Discussion*** 15

### 16 17 18 19 **Microtrap Evaluation** 20

21  
22 This study aimed to simplify an “in-house” thermal desorption (TD) module to perform  
23  
24 environmental air contamination and breath analysis for compounds at sub-ppmv levels  
25  
26 in a reduced time using sample volumes smaller than 1 L. Table 1 shows the list of the  
27  
28 target compounds, chosen as probable smoking biomarkers, used in the microtrap  
29  
30 evaluation process.  
31  
32  
33

34  
35  
36 Commercial instrumentation is based on TD modules containing two different traps. A  
37  
38 preliminary trap is used to retain the compounds of interest. After the adsorption  
39  
40 process, the retained compounds are thermally desorbed but they cannot be directly sent  
41  
42 to a GC column because the bandwidth of the compounds after this thermal desorption  
43  
44 is excessive. Thus, a second cryogenic trap is needed to focalize the compounds. This  
45  
46 cryotrap is also heated and compounds can be sent to a GC column. This process has  
47  
48 two main disadvantages. First, the use of two traps and two heating processes increases  
49  
50 the formation of artifacts and many compounds can be formed during the heating  
51  
52 process, which may lead to the detection of compounds that are not present in the  
53  
54 original sample [20]. Second, the use of a cryotrap increases the complexity of the  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 instrumentation and makes it necessary to use cryogenic gases, which complicates the  
2 development of portable instruments based on this configuration. It has been  
3  
4 demonstrated that it is possible to develop miniaturized TD modules without the use of  
5  
6 cryogenic traps, which allow the direct injection of the retained VOCs in the multibed  
7  
8 trap into a GC column [17,20,21]. When using this simplified configuration, the  
9  
10 formation of artifacts due to thermal oxidation reactions in the surface of the sorbents is  
11  
12 reduced and more reproducible results can be obtained [20].  
13  
14  
15  
16  
17  
18

19 Feng and Mitra [21] showed that the greater the internal diameter of the adsorbent  
20  
21 microtrap, the greater the bandwidth of the desorbed compounds. However, the use of  
22  
23 small-internal diameter traps has the disadvantage that there is an increase in the  
24  
25 pressure drop along the trap, which decreases the sampling flow rate. This represents a  
26  
27 problem when large amount of air samples are to be analyzed. In our case, a trap with a  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

When TD modules without the use of a secondary cryogenic trap are developed, the  
heat transfer through the material of the microtrap support is also important to obtain  
fast heating of carbon-based sorbents. The use of glass tubes for the design of the trap  
simplifies its manipulation and construction. Using glass supports, the heating process  
has to be performed through a heating wire directly connected to the external wall of the  
glass support [22,23]. The use of small, thick layers of glass is recommended to  
facilitate the heat transfer to the carbon sorbents. Unfortunately, the use of glass  
supports results in a moderate heat transfer through the glass wall leading to a slower  
desorption process of the retained compounds in the carbon-based sorbents. As can be  
seen in Figure 2a, the use of a glass support for the trap yields broad peaks at the

1  
2 beginning of the chromatogram making it difficult to separate pairs of compounds 1-2  
3 and 3-4 when the trap is heated at temperatures around 380°C (solid-line plot). Larger  
4 temperatures have to be applied to achieve a faster desorption of the compounds and to  
5 obtain narrow peaks in the chromatograms (broken-line plots). However, when  
6 temperatures above 400°C are applied, carbon-based sorbents are unstable and start to  
7 decompose, resulting in the detection of highly volatile compounds being released from  
8 them (peaks appearing between 30-35 s). The use of metallic supports results in a better  
9 heat transfer to the sorbents, which leads to faster desorption of compounds (Figure 2b).  
10 It is not necessary to apply desorption temperatures above the stability temperature of  
11 the carbon-based sorbents to obtain constant bandwidths for the desorbed compounds.  
12 A desorption temperature around 300°C has been found to be sufficient to obtain  
13 adequate bandwidths [17], which also helps to prevent thermal degradation of the most  
14 reactive VOCs [20].  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33

34 The use of carbon molecular sieve (CMS) sorbents in the multibed trap suffers from the  
35 disadvantage that this bed determines the water uptake of the system, which is a  
36 problem when air samples containing large amounts of water vapour have to be  
37 analyzed, as in the case of breath samples. Unfortunately, CMS are required for the  
38 quantitative adsorption of compounds that are more volatile than pentane and high  
39 polarity low-molecular-weight compounds [24,25]. The evaluation of a multibed trap  
40 without the presence of a CMS showed that some of the main components in breath  
41 samples (e.g., acetone, methanol, ethanol) cannot be adequately adsorbed when only  
42 graphitized carbons are present in the sorbent tube. Given that Carboxen 1000 saturates  
43 at larger volumes of water than other CMS [26], it would seem to be the best choice for  
44 the analysis of breath samples.  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1  
2 Similar trap configurations to the one used in this study have previously been described  
3  
4 using four-bed [17,20] and three-bed configurations [14]. The use of a three-bed trap is  
5  
6 adequate for the analysis of breath samples and a previous study has managed to detect  
7  
8 more than 250 VOCs using this configuration [14]. However, the cost of a GCxGC-  
9  
10 TOF-MS instrument and the complexity of the technique make it unsuitable for routine  
11  
12 analysis in which it is necessary to evaluate a large number of samples in short periods  
13  
14 of time. The use of conventional GC instruments with simple detectors is preferable for  
15  
16 this type of study.  
17  
18  
19  
20  
21  
22  
23

24 GC with FID as a detector is the most common and simplest instrumentation for the  
25  
26 analysis of VOCs. FID is a suitable detector when a synthetic sample with a known  
27  
28 composition is analyzed. However, when a complex sample such as breath is evaluated,  
29  
30 this detector yields many peaks that cannot be correctly assigned to a specific  
31  
32 compound because of the complexity and variability of samples. As can be seen in  
33  
34 Figure 3, there are many compounds detected in a breath sample at ppbv-pptv levels. In  
35  
36 the case illustrated, only acetone, isoprene (the main components in all breath samples),  
37  
38 and hexane (a known environmental air contaminant in this sample) can be adequately  
39  
40 assigned due to the lack of specificity of the detector. Unfortunately, the use of an FID  
41  
42 only allows screening analysis of breath samples, making it impossible to perform an  
43  
44 adequate study of biomarkers using a GC-FID instrument. An MS detector is more  
45  
46 suitable for the performance of efficient biomarker studies in breath samples.  
47  
48  
49  
50  
51  
52  
53  
54  
55

56 Figure 4 shows the chromatograms of the analysis of breath samples from a smoker and  
57  
58 a non-smoker using the MS detector in scan mode with an  $m/z$  interval from 35 to 200.  
59  
60  
61  
62  
63  
64  
65

1 The target compounds were quantified using a target ion and identified using qualifier  
2 ions and retention times, which are shown in Table 1. The large variation in the  
3  
4 concentrations of the compounds in breath samples (acetone at ppmv level and other  
5  
6 target VOCs at ppbv and pptv levels) results in a non-smoker chromatogram in which  
7  
8 many compounds are seen as minor components that can easily be associated with the  
9  
10 noise of the chromatogram when poorly trained staff perform the analysis.  
11  
12  
13  
14  
15  
16  
17  
18

### 19 **Quantitative Analysis**

20 Calibration data was evaluated for all the target compounds selected. Calibration ranges  
21  
22 (5 ng to 150  $\mu\text{g}$ ) were chosen to cover a mass range similar to that expected in breath  
23  
24 samples for healthy non-smokers. Sampling times similar to those used in breath  
25  
26 samples analysis (5-25 minutes) were also used in the calibration analysis. Three  
27  
28 replicates were conducted at each sampling time, and five different mass values were  
29  
30 used for calibrations. All compounds showed linearity in the range studied and the  
31  
32 determination coefficients were above 0.99 in all cases, which can be considered  
33  
34 adequate taking into account the small concentrations of the standards and the use of  
35  
36 adsorption times of up to 25 minutes.  
37  
38  
39  
40  
41  
42  
43  
44  
45

46 The quantification limits obtained for all the compounds were in the 20-40 pptv range  
47  
48 for 775  $\text{cm}^3$  of sample. The experimental detection limits (determined for peaks with a  
49  
50 signal-to-noise ratio of five in the chromatograms) were between 5-10 pptv for the  
51  
52 target compounds. The use of larger sample volumes could lead to reduced  
53  
54 quantification and detection limits. However, this would also increase the analysis time  
55  
56 and require a large amount of sample. The results obtained in the analysis of 31  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 different breath samples from smokers and non-smokers show that the quantification  
2 limits obtained in the analysis of 775 cm<sup>3</sup> of breath sample was enough to determine all  
3 the target compounds evaluated. The analysis of sample volumes below 620 cm<sup>3</sup> does  
4 not allow the identification of one of the target compounds, 2,5-dimethylfuran, in some  
5 samples. In the chosen conditions, 40 minutes is sufficient to perform the complete  
6 analysis (adsorption, thermal desorption, and GC analysis).  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18

### 19 **Environmental Air and Breath Analysis**

20  
21  
22  
23

24 The objective of the device being developed in this study is not to determine all the  
25 compounds present in breath samples but rather to determine the specific target  
26 compounds set out in Table 1.  
27  
28  
29  
30  
31  
32  
33

34 As indicated in the “Trap Evaluation” section, the use of FID detection does not allow  
35 the correct identification of the target VOCs in breath samples. However, the use of this  
36 detector has been useful to observe the effect of environmental air contamination in the  
37 breath of volunteers working in different laboratories. Figure 5 shows the  
38 chromatogram of a volunteer non-smoker working in a chemistry laboratory (Fig. 5a).  
39 The peaks detected at around 3, 4, and 6 minutes are not usually present at these levels  
40 in the breath of healthy volunteers. Posterior analysis of the environmental air of the  
41 laboratory where the volunteer was working (Fig. 5b) showed the presence of these  
42 compounds in the air of the laboratory, which were also breathed by the volunteer  
43 resulting in exogenous contamination of his breath. The analysis of a sample from  
44 another volunteer in the same laboratory after the ingestion of a sweet shows the  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1 presence of the same environmental air contamination as well as the exogenous  
2 contamination from the flavours in the sweet (Fig. 5c).  
3  
4  
5

6  
7 The contamination due to exposure in working environments was evaluated during  
8 working and non-working periods. Figure 6 shows the evolution of the concentration  
9 profiles for four VOCs (acetone, methylene chloride, hexane, and ethyl acetate) in the  
10 environmental air of a chemistry laboratory during a working day. These solvents are  
11 used daily in an adjacent laboratory where basic safety precautions (e.g., working in  
12 safety cabinets with fume hoods) are taken during their manipulation. However, some  
13 poor working habits were also revealed (e.g., laboratory doors were usually open and  
14 VOCs had the possibility to be diffused into common areas and other laboratories). The  
15 routine use of acetone and methylene chloride in this laboratory releases higher contents  
16 of these VOCs into the air during working hours. Profiles obtained for hexane and ethyl  
17 acetate show some undetermined maximum levels at specific time periods as these  
18 solvents were used discontinuously. Note that at the beginning of the working day  
19 levels detected for all compounds are minimum because these solvents were not in use  
20 during the night. A slight decrease in the contamination was also observed just after  
21 13:00 h due to the lunch break. The analysis of the environmental air of the same  
22 laboratory during non-working days was always significantly lower, with levels  
23 detected above detection limits but below quantification (concentrations ranging  
24 between 10-40 pptv).  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52

53 These results are very promising as they demonstrate the possibility of the developed  
54 system performing near real-time air contamination analyses. Other similar systems,  
55 and also common passive sampling systems, only allow the determination of the  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65



1 average contamination over a period of some hours. The developed system permits one  
2 analysis to be made every 40 minutes with detection limits in the low pptv levels  
3  
4 ( $\text{ng}\cdot\text{m}^{-3}$ ). This leads to a more precise knowledge of the concentration profile of  
5  
6 contaminants and may help for preventative measures to be taken faster possibly  
7  
8 avoiding dangerous conditions.  
9

10  
11  
12  
13  
14 The analysis of employees' breath samples working in those laboratories where  
15  
16 contamination was found showed these contaminants to be present in their breath at  
17  
18 higher levels than were detected in volunteers not working near to a laboratory (Figure  
19  
20 7). In the case of methylene chloride, concentrations detected in breath samples of  
21  
22 employees working in a chemistry laboratory ranged from 4 to 103 ppbv ( $n=16$ ,  $sd=26$ ).  
23  
24 Levels detected for those who do not work in or near a chemistry laboratory ranged  
25  
26 from 0.9 to 4.5 ppbv ( $n=15$ ,  $sd=0.8$ ). The large variability found in the first group is due  
27  
28 to the different exposition of the volunteers to the contaminant. 103.1 ppbv were  
29  
30 detected in a person who was directly manipulating this solvent (using a fume hood).  
31  
32 Other employees who did not manipulate this solvent reached values of between 4 and  
33  
34 49 ppbv, depending on their proximity to the source and the time of day. These results  
35  
36 show that the use of conventional safety precautions (conventional fume hoods) when  
37  
38 working with solvents containing VOCs is not sufficient to prevent contamination of the  
39  
40 environmental air. Moreover, the high diffusivity of the most volatile compounds also  
41  
42 results in the contamination of the surrounding areas. This means that more intensive  
43  
44 safety precautions should be taken when these solvents are used.  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54

55  
56 The results obtained in the evaluation of the target compounds chosen as probable  
57  
58 smoking biomarkers (Table 2) fit well with the findings of earlier studies [8-  
59  
60  
61  
62  
63  
64  
65

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

11,13,14,16]. Benzene, toluene, and *m*-,*p*-xylene concentrations detected in smokers (n=14) show a significant difference ( $p<0.001$ , test U of Mann-Whitney) to those levels found in non-smokers (n=17). The evaluation of the data obtained for *o*-xylene does not show significant difference between the two groups ( $p=0.082$ ). The main problem observed with benzene, toluene, and *m*-,*p*-xylene is that all these compounds are also present in non-smokers and their variability in non-smokers is sufficiently large as to make it difficult to determine their smoking status unless samples from smokers are analyzed shortly after the last puff of the cigarette. As can be seen in Table 2, the maximum concentration levels found for these compounds in non-smokers were always higher than the minimum levels found in smokers. Previous studies [9,12,13] found that benzene content in breath decreases rapidly to pre-exposure levels in less than 1 hour after smoking.

In order to distinguish between endogenous substances and exogenous contaminants, corrections for background concentrations of volatile compounds in inhaled air is mandatory. For this reason, the most promising results for probable smoking biomarkers seem to be obtained with 2,5-dimethylfuran. This compound has only been detected in one sample from a non-smoker, who was a passive smoker living in a house with smokers and had more than 3 hours of contact with the smoke of cigarettes every day. This suggests that correction for background concentration may not be needed for this compound and also that it might be a strong indicator of passive smoking.

## Conclusions

1 The simplicity and high sensitivity of the capillary TD module developed and the ability  
2 to determine the presence of VOCs at sub-ppbv levels in both environmental air and  
3  
4 breath samples in a short period of time (less than 40 minutes) introduces new  
5  
6 perspectives in this type of analysis. The possibility of obtaining near real-time  
7  
8 determination of VOCs has allowed us to demonstrate that contamination by highly  
9  
10 volatile compounds takes place very easily. Routine safety precautions for working with  
11  
12 solvents containing highly volatile compounds are not sufficient for a full prevention  
13  
14 against contamination of the environmental air of the working space and surrounding  
15  
16 areas.  
17  
18  
19  
20  
21  
22  
23

24 The analysis of probable smoking biomarkers in breath samples shows that compounds  
25  
26 such as benzene and toluene can be used for the determination of smoking status but  
27  
28 only when samples are analyzed shortly after smoking. The results obtained with  
29  
30 2,5-dimethylfuran seem to be very promising as this compound has only been detected  
31  
32 in active smokers and one passive smoker.  
33  
34  
35  
36  
37  
38

39 The possible use of 2,5-dimethylfuran as a smoking biomarker has also been proposed  
40  
41 in other previous studies [8,14-16]. However, the reduced number of samples evaluated  
42  
43 in these studies makes them insufficient for adequate statistical analysis and,  
44  
45 furthermore, there is the need for a clinical study to test this hypothesis and the utility of  
46  
47 the compound as a useful smoking status biomarker. This study is currently been  
48  
49 performed in our laboratory and the results will be evaluated after a suitable number of  
50  
51 samples have been analyzed.  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

1  
2  
3 *Acknowledgements*  
4  
5

6 M. Alonso acknowledges the Spanish Ministry of Education for the financial support of  
7  
8 her collaboration research (Becas colaboración curso 2007-2008, orden  
9  
10 ECI/2007/2007). This study has been financed by the University of Girona, project #  
11  
12 003798 (ajuts R+D en el marc de l'IdIBGi 2006).  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

**References**

- [1] A. Amman, A. Schmid, S. Scholl-Burgi, S. Telser, H. Hinterhuber, *Spectroscopy Europe* 17 (2005) 18.
- [2] T.H. Risby, S.F. Solga, *Appl. Phys. B* 85 (2006) 421.
- [3] M. Phillips, J., Herrera, S. Krishnan, M. Zain, J., Greenberg, R.N. Cataneo, J. *Chromatogr. B* 729 (1999) 75.
- [4] M. Phillips, *Anal. Biochem.* 247 (1997) 272.
- [5] D.K. Hatsukami, N.L. Benowitz, S. Rennard, C. Oncken, S.S. Hecht, *Nicotine and Tobacco Research* 8 (2006) 599.
- [6] N.L. Benowitz, *Environ. Health Perspect.* 107 (1999) 349.
- [7] S.S. Hecht, S.G. Carmella, S.E. Murphy, S. Akerkar, K.D. Brunneman, D. Hoffmann, *N. Engl. J. Med.* 329 (1993) 1543.
- [8] S.M. Gordon, *J. Chromatogr.* 511 (1990) 291.
- [9] A. Jordan, A. Hansel, R. Holzinger, W. Lindinger, *Int. J. Mass Spectrom. Ion Processes* 148 (1995) L1.
- [10] L. Wallace, T. Buckley, E. Pellizzari, S. Gordon, *Environ. Health Perspect.* 104 suppl. 5 (1996) 861.
- [11] L. Wallace, *Environ. Health Perspect.* 104 suppl. 6 (1996) 1129.
- [12] W.K. Jo, K.W. Pack, *Environ. Res.* 83 (2000) 180.
- [13] S.M. Gordon, L.A. Wallace, M.C. Brinkman, P.J. Callahan, D.V. Kenny, *Environ. Health Perspect.* 110 (2002) 689.
- [14] J.M. Sanchez, R.D. Sacks, *Anal. Chem.* 78 (2006) 3046.
- [15] M. Libardoni, P.T. Stevens, J.H. White, R. Sacks, *J. Chromatogr. B* 842 (2006) 13.
- [16] J.J.B.N. Van Berkel, J.W. Dalling, G.M. Moeller, R.W.L. Godschalk, E. Moonen, E.F.M. Wouters, F.J. Van Schooten, *J. Chromatogr. B* 861 (2008) 861, 101.

- 1 [17] J.M. Sanchez, R.D. Sacks, *Anal. Chem.* 75 (2003) 978.  
2  
3 [18] J. Dewulf, H. Van Langenhove, *J. Chromatogr. A* 843 (1999) 163.  
4  
5 [19] K.A. Cope, M.T. Watson, W.M. Foster, S.S. Sehnert, T.H. Risby, *J. Appl. Physiol.*  
6  
7 96 (2004) 1371.  
8  
9 [20] J.M. Sanchez, R.D. Sacks, *J. Sep. Sci.* 28 (2005) 22.  
10  
11 [21] C. Feng, S. Mitra, *J. Chromatogr. A* 805 (1998) 169.  
12  
13 [22] C. Lu, E. Zellers, *Anal. Chem.* 73 (2001) 3449.  
14  
15 [23] C. Lu, E. Zellers, *Analyst* 127 (2002) 1061.  
16  
17 [24] K. Dettmer, T. Knobloch, W. Engewald, *Fresenius J. Anal. Chem.* 366 (2000) 70.  
18  
19 [25] K. Dettmer, W. Engewald, *Anal. Bioanal. Chem.* 373 (2002) 490.  
20  
21 [26] P. Fastyn, W. Kornacki, M. Kardas, J. Gawlowski, J. Niedzielski, *Analyst* 128  
22  
23 (2003) 198.  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

## Figure Captions

1  
2  
3  
4  
5 **Figure 1.** Diagram of the capillary multi-bed sorption microtrap / GC system  
6 developed. Valves  $V_1$ ,  $V_2$ , and  $V_3$  are used to control the gas flow direction through the  
7 trap tube. Inset shows details of the three-bed sorption trap design (see text for details).  
8  
9

10  
11  
12  
13  
14 **Figure 2.** Chromatograms of a synthetic mixture containing five highly volatile  
15 compounds to evaluate the desorption performance of the microtrap after a rapid heating  
16 pulse is applied to the trap for the desorption of the adsorbed VOCs. (a) glass support,  
17 (b) metallic support. Desorption temperatures: 380°C (solid lines) and 500°C (dashed  
18 line).  
19  
20  
21  
22  
23  
24  
25  
26  
27

28  
29 **Figure 3.** GC-FID chromatogram of a human breath sample collected from a volunteer.  
30 A sample collection time from a 1 L gas sampling bag of 25.0 min at  $31 \text{ cm}^3 \cdot \text{min}^{-1}$  was  
31 used ( $775 \text{ cm}^3$  breath sample).  
32  
33  
34  
35  
36  
37

38  
39 **Figure 4.** GC-MS chromatograms of the analysis of breath samples from a smoker after  
40 1 minute of smoking (a) and a non-smoker (b). Detection in scan mode ( $m/z$  range 35-  
41 200). Peak numbers correspond to the compounds number in Table 1.  
42  
43  
44  
45  
46  
47

48  
49 **Figure 5.** GC-FID chromatograms of: a) breath sample from a non-smoker working in a  
50 laboratory, b) the environmental contamination of this laboratory, and c) breath sample  
51 of a non-smoker working in the same laboratory after eating a flavored sweet.  
52  
53  
54

55  
56 Keynote: A - Acetone; C - Environmental air contamination; F – Sweet flavors.  
57  
58  
59  
60  
61  
62  
63  
64  
65

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60  
61  
62  
63  
64  
65

**Figure 6.** Concentration profiles found for acetone, methylene chloride, hexane, and ethyl acetate in the environmental air of a chemistry laboratory during a working day.

**Figure 7.** Median and quartiles for the methylene chloride levels detected in breath samples from employees working in a chemistry laboratory and those who do not work in or near a chemistry laboratory ( $p < 0.001$ ).



**Table 1.** Target compounds in chromatographic elution order, and their quantifier and qualifier ions. Retention times ( $t_R$ ) and method detection limits (MDL) are indicated for the GC-MS determination and the analysis of 775 cm<sup>3</sup> of sample.

#	Compound	Quantifier ion	Qualifier ions	$t_R$ (min)	MDL	
					(pptv) <sup>a</sup>	(ng·m <sup>-3</sup> )
1	benzene	78	77, 51, 52	3.24	5	16
2	2,5-dimethylfuran	96	95, 81, 53	3.89	5	20
3	toluene	91	92	5.03	5	19
4	<i>m</i> -, <i>p</i> -xylene	91	106, 105, 77	6.89	10	43
5	<i>o</i> -xylene	91	106, 105, 77	7.07	10	43

a) values determined at 25°C and 760 mm Hg.

**Table 2.** Statistics obtained on analysing the target compounds for 31 breath samples. Concentration values expressed in ppbv (25°C and 760 mm Hg).

Non-smokers (n=17)

	benzene	dimethylfuran	toluene	<i>m-, p-xylene</i>	<i>o-xylene</i>
Mean	1.94	0.03	1.23	0.10	0.44
Median	0.73	0.00	0.97	0.07	0.23
Variance	5.26	0.02	1.45	0.01	0.52
Minimum	0.12	bdl*	0.14	0.03	0.03
Maximum	6.88	0.58	3.63	0.39	2.60

\* below detection limit (method detection limit = 5 pptv)

Smokers (n=14)

	benzene	dimethylfuran	toluene	<i>m-, p-xylene</i>	<i>o-xylene</i>
Mean	271.87	25.39	89.27	1.00	1.76
Median	28.09	4.73	12.12	0.28	1.19
Variance	41x10 <sup>4</sup>	42 x10 <sup>2</sup>	37 x10 <sup>2</sup>	5.37	1.99
Minimum	5.22	0.37	1.31	0.07	0.44
Maximum	2292.54	246.72	675.76	8.97	5.15

Figure 1




Figure 2




Figure 3





Figure 4







Figure 6





iric

